

ST Hood Short-type

NEW


Small-caliber Endoscopic Hoods



FUJIFILM

- ▶ Maintains an excellent endoscopic field of view during clinical procedures such as ESD

Application of hoods



Sample image (short-type)



- ▶ To be attached to the tips of endoscopes
- ▶ Prevents the operative field from being disturbed by the mucous membrane
- ▶ Ensures a good operative field during endoscopic procedures such as EMR, ESD, POEM, etc.

Features

- ▶ Small-caliber tip short design allows an optimized view to the operative field
- ▶ 7.0 mm distance to the distal end and an inner diameter of 8.0 mm to support efficient visibility
- ▶ Two drains realize the continuous flow of fluid
- ▶ Sterile and single-use packages
- ▶ No guide ditch necessary; enables an easy insertion of accessories

Article number	Outer diameter	Distance from the endoscope tip	Inner diameter of distal end	Drain	Guide ditch
DH-28GR	11.8 mm	7.0 mm	8.0 mm	2	Without
DH-29CR	13.0 mm	7.0 mm	8.0 mm	2	Without
DH-30CR	14.8 mm	7.0 mm	8.0 mm	2	Without

Generic name: Single-use non-active equipment for natural orifice endoscopic treatment. Package: 5 pcs per pack

ST Hood related products

FlushKnife BT


Covering procedures from marking to hemostasis, this slitter equipped with the water jet system offers high versatility. The ball tip produces good traction during dissection.



Generic Name: Active instrument, endoscopic, high-frequency current procedure, single-use

FlushKnife

The water jet system, which cleans the tip of the knife, keeps the knife sharp throughout the operation. This slitter covers from marking to hemostasis.



FUJIFILM

FUJIFILM Europe GmbH

Heesenstr. 31, 40549 Düsseldorf, Germany
Tel.: +49 211-50 89 0, Fax: +49 211-50 89 8700
www.fujifilm.eu