




**C.R.Kennedy**  
| MEDICAL SOLUTIONS |

# FIRST STEP™

Flexible Endoscope Bedside Cleaning Kit



## Ready-To-Use First Step Kits

Convenience of use increases staff compliance and turnover time. The ready-to-use detergent ensures a properly mixed formula every time. The combined effect is 100% compliance, 100% of the time.

## The Add Water First Step Kit

The Add Water Kit contains concentrated Simple2 multi-tiered enzymatic detergent in a convenient stand-up pouch. Ergonomic Krinkle cleaning pad included with contoured ridges to help remove gross contamination.

## Simple2™

Simple2 effectively dissolves blood, fat, tissue, protein and most other forms of organic material. Simple2 is low-foaming, non-toxic, non-corrosive, latex-free, neutral pH and safe to use on all flexible endoscopes.



Before Water      Water Added

The Add Water Kit maximises storage space and reduces freight costs



500ml  
EP-6W

Order Code	Description	Unit
EP-6W	First Step Bedside Kit – 500ml Stand-Up Pouch with Krinkle Cleaning Pad	120/box



Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
medicalcustomerservice@crkennedy.com.au | crkennedy.com.au