



C.R.Kennedy
| MEDICAL SOLUTIONS |

HELICOTEC UT[®] PLUS

Helicobacter Pylori Paper Urease Test

Detects the urease activity of Helicobacter pylori in gastric mucosal biopsies

Clear, rapid, and accurate results within 5 minutes to 1 hour

Convenience of providing medication on the day of the procedure

Superior quality, high sensitivity and specificity

Does not need refrigeration. Can be stored in the procedure room for convenience

Includes a specimen pick for easy transfer

HelicotecUT[®]Plus is stable for up to 30 months



Easy to Use



Peel the label



Transfer the biopsy specimen onto the paper with the specimen pick enclosed



Re-seal and squeeze



Record the date, time, and patient info. Then read the result

Order Code	Description	Unit
HUP01	Helicotec UT [®] Plus Paper Urease Test	50/Box



Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544
medicalcustomerservice@crkennedy.com.au | crkennedy.com.au