

OASIS[®] SCOPE TRANSPORT

Scope Transport Tray with Center Reservoir

- 

Single-Use and disposable.
- 

Made of a biodegradable plant-based material.
- 

Space efficient, stackable and nestable.

Single-use rigid containment

Bedside cleaning is performed within the Oasis Tray, keeping residual run-off contained

Rigid construction protects delicate scopes, while the reversible lid clearly differentiates clean and soiled scopes

Single-use alternative that protects both patients and staff from contamination and infection

Oasis Trays are latex-free



Order Code	Description	Unit
OT1230	Oasis Biodegradable Tray with centre reservoir and reversible clean/soiled lid.	30/box



Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544
 medicalcustomerservice@crkennedy.com.au | crkennedy.com.au