

CLEANING BRUSHES

Premium Single-use Endoscope Channel Cleaning Devices



Reliable cleaning brushes

Rounded brush head designed to assist with reducing scope damage

Available with or without toothbrush

Complete kits to facilitate cleaning room efficiencies

Order Code	Description	Length	∅ Working Channel	Unit
WILBKA	Double Ended Cleaning Brush, Port Hole Cleaning Brush, and Toothbrush	230cm	2.8 - 4.2mm	50/Box
WILBKQ	Double Ended Cleaning Brush, and Port Hole Cleaning Brush	230cm	2.8 - 4.2mm	50/Box
WILBKI	Single Ended Cleaning Brush, Port Hole Cleaning Brush, Single Ended Cleaning Brush, and Toothbrush	100cm 6cm 30cm	2.0 - 2.2mm 10.0mm 2.8 - 4.2mm	50/Box
WILBKH	Single Ended Cleaning Brush, Port Hole Cleaning Brush, Single Ended Cleaning Brush, and Toothbrush	100cm 6cm 30cm	1.2mm 10.0mm 2.8 - 4.2mm	50/Box



Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544
medicalcustomerservice@crkennedy.com.au | crkennedy.com.au