

# ST HOOD

Fujifilm® ST Hood for EMR, ESD and POEM procedures



Ensures a good operative field during endoscopic procedures such as EMR, ESD, POEM, etc

Two drains enhances the continuous flow of fluid

Small-caliber tip short design allows an optimized view

7.0 mm distance to the distal end to support efficient visibility

Sterile and single-use, individually packaged

Order Code	Hood Outer Diameter	Compatible Fujifilm® Endoscopes	Compatible Olympus® Endoscopes	Unit
DH-28GR	11.8mm	EG-760R, EG-760Z, EC-740T	GIF-HQ190, GIF-1TH190, GIF-H180, PCF-PH190	5/Box
DH-29CR	13.0mm	EG-760CT, EC-760P		5/Box
DH-30CR	14.8mm	EC-760R, EC-760ZP, EC-760S, EC-760Z	CF-H190, CF-HQ190, PCF-H190D, PCF-H180A, GIF-2TH180	5/Box



Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
[medicalcustomerservice@crkennedy.com.au](mailto:medicalcustomerservice@crkennedy.com.au) | [crkennedy.com.au](http://crkennedy.com.au)