





C.R.Kennedy
| MEDICAL SOLUTIONS |

TWIST 'n' CATCH

Fujifilm® Medwork Twist 'n' Catch Stone Extraction Basket



- Maximum expansion force
- Precision control by rotary motion
- Infinite adjustment of basket size
- Proven Twist 'n' Catch basket captures stones quickly and reliably
- Spherical 130° rotation for short and long wire technique
- Manufactured in Germany



Maximum expansion power

Order Code	Description	∅ Basket	Guide Wire	Working Length	∅ Working Channel	Unit
E502285	BAS1-K3-35-28-200 Twist n Catch Basket For short and long wire technique	35mm	.035" ≥ 200cm	200cm	3.2mm	5/Box

ATTENTION: Medwork products of the LIT series are perfectly suited for use with other system products in this series. Please note that when combined with products / instruments from other manufacturers substantial functional limitations may occur. A liability of medwork for such cases is excluded.



Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544
medicalcustomerservice@crkennedy.com.au | crkennedy.com.au