

# JAZZ VALVE BUNDLE

The only disposable valve kit recommended by Fujifilm® Japan for Eluxeo™ 700 Series endoscopes combined with a Water Jet Adapter with Non-Return Valve




Manufactured by Fujifilm® Medwork, Germany  
 Protection against infection and cross contamination  
 Reduce time and effort in endoscope reprocessing  
 Perfect valve operation - every time  
 Disposable endoscope valve kit consisting of air/water valve, suction valve and biopsy valve FROG  
 For Fujifilm® endoscopes (G7 series), compatible with ED-580XT and EI-580BT, and 700 Series endoscopes

Help to eliminate the risk of cross-contamination from backflow into the irrigation connector

Includes backflow valve

Sterile and designed as single use



Order Code	Description	Unit
P3-700	VAL1-G7-50 Air/Water, Suction and Biopsy Single-use Valves for Fujifilm® 700 Series / G7 With BKNRF Auxiliary Water Jet Connector for Fujifilm® Endoscopes	100 Sets/Box



Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
 medicalcustomerservice@crkennedy.com.au | crkennedy.com.au