



C.R.Kennedy
| MEDICAL SOLUTIONS |

INJECTION NEEDLE

Single use injection needle for endoscopy



CE certified Disposable Injection Needle
Exceptionally sharp precision cut needle
Compatible with all endoscope manufacturer brands
Very easy to use with a smooth operation
Stiff outer PTFE tube to prevent kinking



Order Code	Description	Working Length	∅ Working Channel	Unit
WS-2423PN2504	Single Use injection needle 4mm 25G	230cm	2.8mm	10/Box
WS-2423PN2505	Single Use injection needle 5mm 25G	230cm	2.8mm	10/Box
WS-2423PN2304	Single Use injection needle 4mm 23G	230cm	2.8mm	10/Box
WS-2423PN2305	Single Use injection needle 5mm 23G	230cm	2.8mm	10/Box



Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544
medicalcustomerservice@crkennedy.com.au | crkennedy.com.au