

BIOPSY FORCEPS

Designed to maximise sampling and streamline procedures



Maximise your tissue sampling with the high volume cup

Catheter markers assist with quick insertion and withdrawal*

Hinged jaws aid positioning at the sampling site

Responsive closing force in tortuous anatomy



Order Code	Description	Cup Type	∅ Cup	Working Length	∅ Working Channel	Unit
WS-2423BUZ	Disposable Biopsy Forceps	Serrated	2.4mm	230cm	2.8mm	10/Box
WS-2423BNZ	Disposable Biopsy Forceps with needle	Serrated	2.4mm	230cm	2.8mm	10/Box
WS-2423BTZ	Disposable Biopsy Forceps	Oval	2.4mm	230cm	2.8mm	10/Box
WS-2423BWZ	Disposable Biopsy Forceps with needle	Oval	2.4mm	230cm	2.8mm	10/Box
WS-3023BTZ	Disposable Jumbo Biopsy Forceps	Oval	3.0mm	230cm	3.2mm	10/Box
WS-1816BT	Disposable Mini Biopsy Forceps	Oval	1.8mm	160cm	2.0mm	10/Box
WS-1810BT	Disposable Bronchoscope Biopsy Forceps	Oval	1.8mm	100cm	2.0mm	10/Box

* Available on 2.4mm forceps



Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544
medicalcustomerservice@crkennedy.com.au | crkennedy.com.au