


HELICLIP

Haemostatic clip for endoscopic marking and haemostasis



Repeatable opening and closing
 180° opening and 360° stepless bidirectional rotation
 PTFE sheath coating for optimal sliding in the scope
 Atraumatic tip protects the scope from damage
 Short clip stem length for multiple clip placements

Order Code	Description	Opening Span	Opening Angle	Rotation	Working Length	∅ Working Channel	Unit
HC-23-230-12	Heliclip Rotatable	12mm	180°	360°	230cm	2.8mm	5/Box
HC-23-230-16	Heliclip Rotatable	16mm	180°	360°	230cm	2.8mm	5/Box



Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544
 medicalcustomerservice@crkennedy.com.au | crkennedy.com.au