

C.R. Kennedy  
MEDICAL DIVISION

# PRODUCT GUIDE


## Endoscopy Accessories




Biopsy Forceps



Procedural Accessories



ERCP



Cleaning, Sterilisation & Disinfection



Haemostasis



Polypectomy, EMR & ESD



Distal Caps & Hoods

2022

Procedural Accessories	<b>3</b>
Distal Caps & Hoods	<b>20</b>
Polypectomy, EMR & ESD	<b>23</b>
Haemostasis	<b>26</b>
Cleaning, Sterilisation & Disinfection	<b>32</b>
ERCP	<b>39</b>
Biopsy Forceps	<b>42</b>
How to Order	<b>44</b>



C.R. Kennedy  
MEDICAL DIVISION

# Procedural Accessories



## Single-use Valve Sets for FujiFilm® Endoscopes

- German precision designed for a perfect fit every time
- Compatible with FujiFilm (R) 700, ED-580XT and EI-580BT series endoscopes
- Sterile single use valves to reduce cross-contamination risk and eliminate reprocessing time



**Biopsy Valve**



**Aux Water Jet Connector**



**Suction Valve**



**Air/Water Valve**



**Cleaning Adapter**

Order Code	Description	Unit
ADA-02-JC-FJ	Aux Water Jet Connector for FujiFilm® 700 series endoscopes	100 / BOX
AWV-KIT-07-FJ	Cleaning Valve for FujiFilm® 700 series endoscopes	100 / BOX
AWV-KIT-01-FJ	Disposable 3 piece Endoscope Valve Set for FujiFilm® 700 series endoscopes with Biopsy Valve, Air-Water Valve, Suction Valve	25 / BOX
AWV-KIT-05-FJ	Disposable 4 piece Endoscope Valve Set for FujiFilm® 700 series endoscopes with Biopsy Valve, Air-Water Valve, Suction Valve, Aux Water Jet Connector	25 / BOX
AWV-KIT-08-FJ	Disposable 5 piece Endoscope Valve Set for FujiFilm® 700 series endoscopes with Biopsy Valve, Air-Water Valve, Suction Valve, Cleaning Valve, Aux Water Jet Con.	25 / BOX

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
[medicalcustomerservice@crkennedy.com.au](mailto:medicalcustomerservice@crkennedy.com.au) | [crkennedy.com.au](http://crkennedy.com.au)

## Single-use Valve Sets for Olympus® Endoscopes

- German precision designed for a perfect fit every time
- Compatible with Olympus 180,190 and EVIS X1 series endoscopes
- Sterile single use valves to reduce cross-contamination risk and eliminate reprocessing time



**Biopsy Valve**



**Aux Water Jet Connector**



**Suction Valve**



**Air/Water Valve**



**Cleaning Adapter**

Order Code	Description	Unit
ADA-JC-OUS	Aux Water Jet Connector for Olympus® series endoscopes	100 / BOX
AWV-KIT-07-OLY	Cleaning Valve for Olympus® series endoscopes	100 / BOX
AWV-KIT-02-OLY	Disposable 3 piece Endoscope Valve Set for Olympus® series endoscopes with Biopsy Valve, Air-Water Valve, Suction Valve	25 / BOX
AWV-KIT-01-OLY	Disposable 4 piece Endoscope Valve Set for Olympus® series endoscopes with Biopsy Valve, Air-Water Valve, Suction Valve, Aux Water Jet Connector	25 / BOX
AWV-KIT-05-OLY	Disposable 5 piece Endoscope Valve Set for Olympus® series endoscopes with Biopsy Valve, Air-Water Valve, Suction Valve, Cleaning Valve, Aux Water Jet Con.	25 / BOX

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
[medicalcustomerservice@crkennedy.com.au](mailto:medicalcustomerservice@crkennedy.com.au) | [crkennedy.com.au](http://crkennedy.com.au)


## Single-use valve kit recommended by FujiFilm® Japan for Eluxeo™ 700 Series Endoscopes

- Manufactured by FujiFilm® Medwork, Germany
- Sterile single use valves to reduce cross-contamination risk and eliminate reprocessing time
- Disposable endoscope valve kit consisting of air/water valve, suction valve and biopsy valve FROG
- Compatible with FujiFilm 700 series endoscopes, ED-580XT and EI-580BT



Order Code	Description	Unit
E502562	3 Piece Air/Water, Suction and Biopsy Valve Set for FujiFilm® 700 Series	50 Sets/Box

### Kit Components



Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
[medicalcustomerservice@crkennedy.com.au](mailto:medicalcustomerservice@crkennedy.com.au) | [crkennedy.com.au](http://crkennedy.com.au)

# CLEANING VALVE

## Air/Water Valve for FujiFilm 700 Series Endoscopes

- Manufactured by Body Products RELAX, GmbH
- German precision for a perfect fit every time
- Sterile single use valves to reduce cross-contamination risk and eliminate reprocessing time



Order Code	Description	Unit
AWV-KIT-07-FJ	Cleaning Air/Water Valve (Single Use for FujiFilm 700 Series)	100 Sets/ Box



Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
medicalcustomerservice@crkennedy.com.au | crkennedy.com.au

# CLEANING VALVE

## Air/Water Valve for Olympus Endoscopes

- Manufactured by Body Products RELAX, GmbH
- Compatible with Olympus endoscopes
- German precision designed for a perfect fit every time
- Sterile single use valves to reduce cross-contamination risk and eliminate reprocessing time



Order Code	Description	Unit
AWV-KIT-07-OLY	Cleaning Air/Water Valve (for Olympus Endoscopes)	100 Sets/Box



Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
medicalcustomerservice@crkennedy.com.au | crkennedy.com.au





## Suction Cleaning Adaptor


- Suitable for use with all major brands of endoscopes
- Sterile single use valves to reduce cross-contamination risk and eliminate reprocessing time
- For use in aspirating water and enzymatic fluids from the distal end of the endoscope through the biopsy port and instrument channel



Order Code	Description	Optional Details	Unit
OLY-CT-001	Universal Cleaning Adaptor Tube for Biopsy Port	Compatible with all endoscopes	50 / BOX

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
medicalcustomerservice@crkennedy.com.au | crkennedy.com.au



### For sterile water bottles with CO2 connection

- Manufactured by Body Products RELAX, GmbH
- Sterile disposable tubing for up to 24-hours use
- With hose set and screw cap for sterile Baxter or Braun water bottles

Order Code	Description	Optional Details	Unit
WB-FJ3-WC2	Ecoflush Cap FujiFilm 500/600	Braun / Fresenius	10 / BOX
WBB-FJ3-WC2	Ecoflush Cap FujiFilm 500/600	Baxter	10 / BOX
WB-FJ6-WC2	Ecoflush Cap FujiFilm 700 series	Braun / Fresenius	10 / BOX
WBB-FJ6-WC2	Ecoflush Cap FujiFilm 700 series	Baxter	10 / BOX
WB-OL3-WC2	Ecoflush Cap Olympus	Braun / Fresenius	10 / BOX
WBB-OL3-WC2	Ecoflush Cap Olympus	Baxter	10 / BOX
WB-PX6-WC2	Ecoflush Cap Pentax	Braun / Fresenius	10 / BOX
WBB-PX6-WC2	Ecoflush Cap Pentax	Baxter	10 / BOX

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
medicalcustomerservice@crkennedy.com.au | crkennedy.com.au

## With CO2 Connection Tube

- Manufactured by Body Products RELAX GmbH
- For sterile water bottles where CO2 delivery is required
- Sterile disposable tubing for up to 24-hours use
- Compatible with Baxter and Braun/Fresenius water bottles
- With 1.1m connecting CO2 tube with luer lock connection
- Adaptor for FujiFilm GW-100 Insufflation pump included
- Adaptor for Olympus UCR Insufflator sold separately




Order Code	Description	Optional Details	Unit
WB-FJ4-WC2	Ecoflush Cap with CO2 connection tube FujiFilm 500/600	Braun / Fresenius	15 / BOX
WBB-FJ4-WC2	Ecoflush Cap with CO2 connection tube FujiFilm 500/600	Baxter	15 / BOX
WB-FJ5-WC2	Ecoflush Cap with CO2 connection tube FujiFilm 700	Braun / Fresenius	15 / BOX
WBB-FJ5-WC2	Ecoflush Cap with CO2 connection tube FujiFilm 700	Baxter	15 / BOX
WB-OL4-WC2	Ecoflush Cap with CO2 connection tube Olympus	Braun / Fresenius	15 / BOX
WBB-OL4-WC2	Ecoflush Cap with CO2 connection tube Olympus	Baxter	15 / BOX
WB-PX5-WC2	Ecoflush Cap with CO2 connection tube Pentax	Braun / Fresenius	15 / BOX
WBB-PX5-WC2	Ecoflush Cap with CO2 connection tube Pentax	Baxter	15 / BOX

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
[medicalcustomerservice@crkennedy.com.au](mailto:medicalcustomerservice@crkennedy.com.au) | [crkennedy.com.au](http://crkennedy.com.au)

## ECOFLUSH Cap with CO2 Connection Tube Plus Inline Filter

- Manufactured by Body Products RELAX GmbH
- For sterile water bottles where CO2 delivery is required
- Sterile disposable tubing for up to 24-hours use
- Compatible with Baxter and Braun/Fresenius water bottles
- With separate 1.1m connecting CO2 tube with luer lock connection
- Adaptor for FujiFilm GW-100 Insufflation pump included
- Adaptor for Olympus UCR Insufflator sold separately
- Optional inline filter for patient and insufflation unit protection



Order Code	Description	Optional Details	Unit
WB-FJ3-TU-210	Ecoflush Cap with CO2 connection tube plus inline filter FujiFilm 500/600	Braun / Fresenius	8 / BOX
WBB-FJ3-TU-210	Ecoflush Cap with CO2 connection tube plus inline filter FujiFilm 500/600	Baxter	8 / BOX
WB-FJ6-TU-210	Ecoflush Cap with CO2 connection tube plus inline filter FujiFilm 700	Braun / Fresenius	8 / BOX
WBB-FJ6-TU-210	Ecoflush Cap with CO2 connection tube plus inline filter FujiFilm 700	Baxter	8 / BOX
WB-OL3-TU-210	Ecoflush Cap with CO2 connection tube plus inline filter Olympus	Braun / Fresenius	8 / BOX
WBB-OL3-TU-210	Ecoflush Cap with CO2 connection tube plus inline filter Olympus	Baxter	8 / BOX
WB-PX6-TU-210	Ecoflush Cap with CO2 connection tube plus inline filter Pentax	Braun / Fresenius	8 / BOX
WBB-PX6-TU-210	Ecoflush Cap with CO2 connection tube plus inline filter Pentax	Baxter	8 / BOX

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544

medicalcustomerservice@crkennedy.com.au | crkennedy.com.au

## ECOFLUSH Hybrid CO2 Cap

- Manufactured by Body Products RELAX GmbH
- Unique 3 in 1 cap design
- For sterile water bottles where CO2 delivery and irrigation delivery is required
- Sterile disposable tubing for up to 24-hours use
- Compatible with Baxter and Braun/Fresenius water bottles
- Includes CO2 connection tube and twin irrigation tube from a single cap
- Adaptor for FujiFilm GW-100 Insufflation pump included
- Adaptor for Olympus UCR Insufflator sold separately



Order Code	Description	Optional Details	Unit
WB-HYB-FJ4	HybriX EcoFlushCap Waterbottle Cap with Co2 tubingsystem and Twin Irrigation tube for FujiFilm 500/600 series endoscopes	Braun / Fresenius	10 / BOX
WBB-HYB-FJ4	HybriX EcoFlushCap Waterbottle Cap with Co2 tubingsystem and Twin Irrigation tube for FujiFilm 500/600 series endoscopes	Baxter	10 / BOX
WB-HYB-FJ5	HybriX EcoFlushCap Waterbottle Cap with Co2 tubingsystem and Twin Irrigation tube for FujiFilm 700 series endoscopes	Braun / Fresenius	10 / BOX
WBB-HYB-FJ5	HybriX EcoFlushCap Waterbottle Cap with Co2 tubingsystem and Twin Irrigation tube for FujiFilm 700 series endoscopes	Baxter	10 / BOX
WB-HYB-OL4	HybriX EcoFlushCap Waterbottle Cap with Co2 tubingsystem and Twin Irrigation tube for OLYMPUS endoscopes	Braun / Fresenius	10 / BOX
WBB-HYB-OL4	HybriX EcoFlushCap Waterbottle Cap with Co2 tubingsystem and Twin Irrigation tube for OLYMPUS endoscopes	Baxter	10 / BOX
WB-HYB-PX5	HybriX EcoFlushCap Waterbottle Cap with Co2 tubingsystem and Twin Irrigation tube for all PENTAX endoscopes series	Braun / Fresenius	10 / BOX
WBB-HYB-PX5	HybriX EcoFlushCap Waterbottle Cap with Co2 tubingsystem and Twin Irrigation tube for all PENTAX endoscopes series	Baxter	10 / BOX

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544

medicalcustomerservice@crkennedy.com.au | crkennedy.com.au

## Cost-effective disposable irrigation tubing



- Sterile disposable tubing for up to 24-hours use
- Twin tube design for greater compatibility to all common flushing pumps
- Superior water flow rate
- For use with FujiFilm®, Olympus® and Pentax® endoscopes (with compatible auxilliary water jet adaptors)

Order Code	Description	Bottle Compatibility	Package
IT-OUS-210	Ecoflush 24hr Twin Irrigation Tubing	Braun / Fresenius	10 / BOX
ITB-OUS-210-BX	Ecoflush 24hr Twin Irrigation Tubing	Baxter	10 / BOX

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
medicalcustomerservice@crkennedy.com.au | crkennedy.com.au



## Disposable CO<sub>2</sub> insufflation tubing

- Sterile disposable tubing for up to 24-hours use
- Connects to any CO<sub>2</sub> water bottle cap utilising a lure connector
- Compatible with any CO<sub>2</sub> insufflator with a lure connector or lure adaptor
- Filter designed to offer both machine and patient protection



Order Code	Description	Unit
CO2-TU-210	Ecoflush CO <sub>2</sub> tubing with inline filter	10 / BOX

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
medicalcustomerservice@crkennedy.com.au | crkennedy.com.au

# BITE BLOCK

Oxyshield™ Endoscopic Bite Block



- A simple, safe and highly effective method of quickly delivering oxygen to the patient
- Two nasal and two oral outlets, with side inlet
- Retaining slopes designed to pull the teeth towards the front of the block to help prevent the block from exiting the mouth during the procedure
- Finger friction pad where the nurse may place a finger during the procedure
- 60fr large lumen for easy passage of scopes and dilator bougies
- High density structure to maximise scope protection

Order Code	Description	Unit
2070231K	Oxyshield™ Bite Block with Strap	50/Box
2070231	Oxyshield™ Bite Block	50/Box
2070232	Strap for Bite Block	50/Box




Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
medicalcustomerservice@crkennedy.com.au | crkennedy.com.au



## Endoscope Distal End Protectors

- Expandable mesh protects the tip and lens of the instrument from damage
- Semi-rigid weave allows moisture to escape, preventing potential bacterial growth that may occur with other products
- Fits securely and gentle mesh design will not scratch any part of the instrument or lens



Order Code	Description	Unit
BKPR	Scope Tip Protector - Small Red 2-10mm	100 / BOX
BKPB	Scope Tip Protector - Large Blue 7-15mm	100 / BOX

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
medicalcustomerservice@crkennedy.com.au | crkennedy.com.au



## Endoscopic Marker for marking polyps and other lesions in the gastrointestinal tract

- Immediately ready to use on opening
- Pre-Loaded with 5ml of Black Dye
- Patient tracking labels included in each tube


- Non-toxic formulation containing no India ink
- Sterile, non-pyrogenic, purified carbon formulation
- Black Eye™ has been successfully used in over 50 countries worldwide



Order Code	Description	Unit
TS205	Black Eye™ Endoscopic Marker	10/Box

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
medicalcustomerservice@crkennedy.com.au | crkennedy.com.au



## Helicobacter Pylori Paper Urease Test


- Detects the urease activity of Helicobacter pylori in gastric mucosal biopsies
- Clear, rapid, and accurate results within 5 minutes to 1 hour
- Convenience of providing medication on the day of the procedure
- Superior quality, high sensitivity and specificity
- Does not need refrigeration. Can be stored in the procedure room for convenience
- Includes a specimen pick for easy transfer
- HelicotecUT<sup>®</sup>Plus is stable for up to 30 months



Peel the label and Transfer the biopsy specimen onto the paper with the specimen pick enclosed



Re-seal and squeeze




Record the date, time, and patient info. Then read the result

Order Code	Description	Unit
HUP01	Helicotec UT <sup>®</sup> Plus Paper Urease Test	50/Box

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544

medicalcustomerservice@crkennedy.com.au | crkennedy.com.au



C.R. Kennedy  
MEDICAL DIVISION

# Distal Caps & Hoods



## CHANNEL and 4x SIDE HOLES

- Improve your endoscopic VISION with rapid fluid discharge
- Compliant material for easy attachment with reduced scope and tissue trauma
- Reduced halation compared to clear caps

### How to attach:

Attach the VISION cap to the scope in the direction shown by the arrow. Avoid backwards attachment. The attachment end is softer, thinner and has further distance to the start of the channels



Order Code	Type Name	Scope OD	Compatible FujiFilm® Endoscopes	Compatible Olympus® Endoscopes	Unit
TA02013	SHS	8.6-9.8mm	EG-760R	PCF-PH190	5/Box
TA02014	SHM	9.8-10.6mm	EG-760Z, EG-760CT, EC-740T	GIF-HQ190, GIF-1TH190	5/Box
TA02015	SHL	10.6-12.6mm	EC-760R, EC-760ZP, EC-760P	PCF-H190D, GIF-2TH180	5/Box
TA02016	SHLL	12.6-14.0mm	EC-760S, EC-760Z	CF-H190, CF-HQ190	5/Box
TA02019	SHM-	9.8-10.6mm	EG-760Z, EG-760CT, EC-740T	GIF-HQ190, GIF-1TH190	5/Box



Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
[medicalcustomerservice@crkennedy.com.au](mailto:medicalcustomerservice@crkennedy.com.au) | [crkennedy.com.au](http://crkennedy.com.au)

## FujiFilm® ST Hood for EMR, ESD and POEM Procedures


- Ensures a good operative field during endoscopic procedures such as EMR, ESD, POEM, etc.
- Two drains enhance the continuous flow of fluid
- Small-Calibre Tip, short design allows an optimized view
- 7.0 mm distance to the distal end to support efficient visibility
- Sterile and single-use, individually packaged



Order Code	Hood Outer Diameter	Compatible FujiFilm® Endoscopes	Compatible Olympus® Endoscopes	Unit
DH-28GR	11.8mm	EG-760R, EG-760Z, EC-740T	GIF-HQ190, GIF-1TH190, GIF-H180, PCF-PH190	5/Box
DH-29CR	13.0mm	EG-760CT, EC-760P		5/Box
DH-30CR	14.8mm	EC-760R, EC-760ZP, EC-760S, EC-760Z	CF-H190, CF-HQ190, PCF-H190D, PCF-H180A, GIF-2TH180	5/Box

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
[medicalcustomerservice@crkennedy.com.au](mailto:medicalcustomerservice@crkennedy.com.au) | [crkennedy.com.au](http://crkennedy.com.au)



C.R. Kennedy  
MEDICAL DIVISION

Polypectomy, EMR & ESD



## FujiFilm® Clutch Cutter Therapeutic Forceps for ESD

- 3 in 1 ESD tool for efficient and safe therapeutic procedures
- Insulated outer edge to avoid damaging tissue
- Rotatable distal jaws enable a precise lesion approach
- Toothed jaws grip mucosa membrane securely and efficiently

### Incision

Grip the mucosa by using the toothed jaws, liftup and start with the incision.

### Dissection

Grip and lift up the submucosa, look to confirm, and dissect.

### Coagulation

An efficient coagulation of bleeding spots is possible by using the high frequency jaws.



Order Code	Description	Jaw Length	Working Length	Ø Working Channel	Unit
DP2618DT-35	FujiFilm Single Use Clutch Cutter Scissors	3.5mm	180cm	2.8mm	Individual
DP2618DT-50	FujiFilm Single Use Clutch Cutter Scissors	5.0mm	180cm	2.8mm	Individual

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
[medicalcustomerservice@crkennedy.com.au](mailto:medicalcustomerservice@crkennedy.com.au) | [crkennedy.com.au](http://crkennedy.com.au)



## FujiFilm® Ball-Tip FlushKnife Diathermic Slitter for ESD

### Enhanced suction capability

The ball-tip slim type has a thinner sheath, making the space between the accessory channel and FlushKnife enlarged, which enables increased suction performance while stabilising the sheath

### Excellent manoeuvrability and stability

FlushKnife is stabilised by the minimised gap within the forceps channel

### Improved durability and insertion ability

Reduced resistance during scope insertions and improved durability (less kinking) is achieved by enhanced flexibility and a thinner sheath with the ball-tip slim type knife



### Enhanced Suction and Stability




The ball-tip slim type sheath tip shape is modified compared to a conventional type. The thinner sheath provides increased suction capability, while the shape of the sheath stabilises the axis.

Order Code	Description	Knife Length	Working Length	Ø Working Channel	Unit
DK2620JIB15	FlushKnife BTS Ball Tip Slim Type	1.5mm	200cm	2.8mm	Individual
DK2620JIB20	FlushKnife BTS Ball Tip Slim Type	2.0mm	200cm	2.8mm	Individual
DK2620JIB25	FlushKnife BTS Ball Tip Slim Type	2.5mm	200cm	2.8mm	Individual
DK2620JIB30	FlushKnife BTS Ball Tip Slim Type	3.0mm	200cm	2.8mm	Individual
DK2620JIB15	FlushKnife BT Ball Tip	1.5mm	180cm	2.8mm	Individual
DK2620JIB20	FlushKnife BT Ball Tip	2.0mm	180cm	2.8mm	Individual
DK2620JIB25	FlushKnife BT Ball Tip	2.5mm	180cm	2.8mm	Individual
DK2623JIB15	FlushKnife BT Ball Tip Colonoscope Length	1.5mm	230cm	2.8mm	Individual
DK2623JIB20	FlushKnife BT Ball Tip Colonoscope Length	2.0mm	230cm	2.8mm	Individual

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544

medicalcustomerservice@crkennedy.com.au | crkennedy.com.au




## Designed for easy foreign body and polyp retrieval

- Flexible net design to allow capture and prevent dislodgement
- Available in a wide range of sizes for various clinical needs
- Cost effective to enable a net to be used every time it is needed

Order Code	Description	Ø Net	Working Length	Ø Working Channel	Unit
WS-2423GW25	Retrieval Net	25mm	230cm	2.8mm	10/Box
WS-2423GW30	Retrieval Net	30mm	230cm	2.8mm	10/Box
WS-2423GW40	Retrieval Net	40mm	230cm	2.8mm	10/Box

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
medicalcustomerservice@crkennedy.com.au | crkennedy.com.au



C.R. Kennedy  
MEDICAL DIVISION

# Haemostasis



## Haemostatic clip for endoscopic marking and haemostasis

- Repeatable opening and closing for accurate placement
- Wide opening and fluent 360 degree rotation
- PTFE sheath coating for optimal sliding in the scope
- Atraumatic tip protects the scope from damage
- Short clip stem length for improved manoeuvring and positioning of additional clips
- Coloured handles to distinguish clips sizes, for various clinical needs



Order Code	Description	Opening Span	Claw Angle	Rotation	Working Length	Ø Working Channel	Unit
HC-23-230-12	Heliclip Rotatable	12mm	180°	360°	230cm	2.8mm	10/Box
HC-23-230-16	Heliclip Rotatable	16mm	180°	360°	230cm	2.8mm	10/Box

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
[medicalcustomerservice@crkennedy.com.au](mailto:medicalcustomerservice@crkennedy.com.au) | [crkennedy.com.au](http://crkennedy.com.au)

## AGS Medtech Rotatable Endoscopic Haemostatic Clip

- Repeatable opening and closing for accurate placement
- 360° bidirectional synchronous rotation control
- Excellent jaw closure strength for effective retention
- Simple handle control allows the operator to control the jaw opening span between 11mm and 16mm
- Short clip stem length for improved maneuvering and positioning of additional clips
- Simple sheath design to allow controlled and seamless rotation




Order Code	Description	Opening Span	Claw Angle	Rotation	Working Length	Ø Working Channel	Unit
AG-51043-2300-135-16	Haemostasis Clip AGS	16mm	135°	360°	230cm	2.8mm	10/Box

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
[medicalcustomerservice@crkennedy.com.au](mailto:medicalcustomerservice@crkennedy.com.au) | [crkennedy.com.au](http://crkennedy.com.au)

## Single use injection needle for endoscopy

- Exceptionally sharp precision cut needle and kink-resistant
- Compatible with all endoscope manufacturer brands
- Stiff outer PTFE tube to prevent kinking
- Smooth functioning design even in tortuous applications



Order Code	Description	Gauge	Needle Length	Working Length	Ø Working Channel	Unit
WS-2423PN2504	Single Use injection needle	25G	4mm	230cm	2.8mm	10 / BOX
WS-2423PN2505	Single Use injection needle	25G	5mm	230cm	2.8mm	10 / BOX
WS-2423PN2304	Single Use injection needle	23G	4mm	230cm	2.8mm	10 / BOX
WS-2423PN2305	Single Use injection needle	23G	5mm	230cm	2.8mm	10 / BOX

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
medicalcustomerservice@crkennedy.com.au | crkennedy.com.au

## Single use injection needle for endoscopy


- Manufactured by FujiFilm® Medwork, Germany
- Exceptionally sharp precision cut needle kink-resistant
- Compatible with all endoscope manufacturer brands
- Easy handling due to single hand operation



Order Code	Description	Gauge	Needle Length	Working Length	Ø Working Channel	Unit
INJ1-A2-05-4-23-230	Icus2 Injection Needle 5mm	25G	4mm	230cm	2.8mm	10 / BOX
INJ1-A2-07-3-23-230	Icus2 Injection Needle 7mm	23G	3mm	230cm	2.8mm	10 / BOX
INJ1-A2-07-4-23-230	Icus2 Injection Needle 7mm	23G	4mm	230cm	2.8mm	10 / BOX
INJ1-A2-07-5-23-230	Icus2 Injection Needle 7mm	23G	5mm	230cm	2.8mm	10 / BOX
INJ1-A2-09-5-23-230	Icus2 Injection Needle 9mm	21G	5mm	230cm	2.8mm	10 / BOX
INJ1-A3-05-5-18-240	Manta* Injection Needle 5mm	25G	5mm	240cm	2.8mm	10 / BOX
INJ1-A3-07-5-18-240	Manta* Injection Needle 7mm	23G	5mm	240cm	2.8mm	10 / BOX

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
[medicalcustomerservice@crkennedy.com.au](mailto:medicalcustomerservice@crkennedy.com.au) | [crkennedy.com.au](http://crkennedy.com.au)



C.R. Kennedy  
MEDICAL DIVISION

# Cleaning, Sterilisation & Disinfection





## Premium Endoscope Channel Cleaning Devices

Scientifically proven improvement over existing brush technology



- Unique patented five wiper elements design
- Complete circumferential seal in the lumen tube
- Removes at least 96% of the debris 100% of the time<sup>1</sup>
- A uniform cleaning action = Reduction in deviations in cleaning protocols<sup>2</sup>
- Improve the efficiency and efficacy of the manual cleaning process

Order Code	Description	Colour Code	Length	Ø Working Channel	Unit
2070224	Micro Pull Thru Single Patient Use Pack	●	220cm	1.0 - 1.2mm	60/Box
2070218	Mini Pull Thru Refill Pack	●	220cm	1.4 - 2.6mm	180/Box
2070217	Mini Pull Thru Starter Kit includes Mini Pull Thru, Storage Cassette & Wall Mount	●	220cm	1.4 - 2.6mm	60/Box
2070215	Mini Pull Thru Single Patient Use Pack	●	220cm	1.4 - 2.6mm	60/Box
2070213	Pull Thru Refill Pack	●	220cm	2.8 - 5.0mm	240/Box
2070210	Pull Thru Starter Kit includes Pull Thru, Storage Cassette & Wall Mount	●	220cm	2.8 - 5.0mm	60/Box
2070220	Pull Thru Single Patient Use Pack	●	220cm	2.8 - 5.0mm	60/Box
2070223	Enteroscopy Pull Thru Single Patient Use Pack	●	250cm	2.8 - 5.0mm	60/Box
2070216	Stubby Brush Valve / Control Head Brush (6x30packs)	●	6cm	1.2 - 2.8mm	180/Box
2070211	Push Thru Brush	●	20cm	1.5 - 5.0mm	60/Box

[1] Timothy S Charlton. A comparison of the efficacy of lumen - cleaning devices for flexible gastrointestinal endoscopes. Aust Infect Control 2007; 12(3): 81-90

[2] "...minor deviations from cleaning protocols resulted in persistent microbiological contamination after disinfection" Deva et al. Detection of persistent vegetative bacteria and amplified viral nucleic acid from in-use testing of gastrointestinal endoscopes. J Hosp Inf 1998;39:1 49-157

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544

medicalcustomerservice@crkennedy.com.au | crkennedy.com.au

## Premium Single-use Endoscope Channel Cleaning Devices

- Reliable cleaning brushes
- Rounded brush head designed to assist with reducing scope damage
- Complete kits to facilitate cleaning room efficiencies
- Available with or without toothbrush

Order Code	Description	Unit	Unit	Unit
WILBKA	Double Ended Cleaning Brush, Port Hole Cleaning Brush, and Toothbrush	230cm	2.8 - 4.2mm	50/Box
WILBKQ	Double Ended Cleaning Brush, and Port Hole Cleaning Brush	230cm	2.8 - 4.2mm	50/Box
WILBKI	Single Ended Cleaning Brush, Port Hole Cleaning Brush, Single Ended Cleaning Brush, and Toothbrush	100cm	2.0 - 2.2mm	50/Box
		6cm	10.0mm	
		30cm	2.8 - 4.2mm	
WILBKH	Single Ended Cleaning Brush, Port Hole Cleaning Brush, Single Ended Cleaning Brush, and Toothbrush	100cm	1.2mm	50/Box
		6cm	10.0mm	
		30cm	2.8 - 4.2mm	

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544

medicalcustomerservice@crkennedy.com.au | crkennedy.com.au

## Flexible Endoscope Bedside Cleaning Kit

### Ready-To-Use

Convenience of use improves staff compliance and turnaround time. The ready-to-use detergent ensures a properly mixed formula every time. The combined effect is 100% compliance, 100% of the time

### Add Water

The Add Water Kit contains concentrated Simple2 multi-tiered enzymatic detergent in a convenient stand-up pouch. Ergonomic Krinkle cleaning pad included with contoured ridges to help remove gross contamination

### Simple2™

Simple2 effectively dissolves blood, fat, tissue, protein and most other forms of organic material. Simple2 is low-foaming, non-toxic, non-corrosive, latex-free, neutral pH and safe to use on all flexible endoscopes



The Add Water Kit maximises storage space and reduces freight costs

Order Code	Description	Unit
EP-6W	First Step Bedside Clean Kit – 500ml Stand-Up Pouch with Krinkle Cleaning Pad	120/box

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544

medicalcustomerservice@crkennedy.com.au | crkennedy.com.au

## CAPTURE



Uses split microfibre technology to capture spores and other viable pathogens at a microscopic level. Can help to capture and remove superbug pathogens like C.Diff, VRE, CRE and Norovirus

## REMOVE

Remove gross contamination and provide a more potent disinfection of pathogens. Disposable Wipes can be used with any brand of surface disinfectant and ensure maximum disinfectant contact time

## DISPOSE

Single-use for disposal immediately after use. By capturing and disposing of spores and other viable pathogens, Clear Path wipes prevent the cross contamination of surface areas



Split microfibre 1/100th the thickness of a human hair

Convenient waterproof dispenser package



Order Code	Description	Unit
W5412-30	Clear Path disposable lint-free low-particulate absorbent microfibre wipes 30 wipes per waterproof dispenser package, 20 packages per box	600 wipes/ box

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544

medicalcustomerservice@crkennedy.com.au | crkennedy.com.au

## SCOPE TRANSPORT BAG

- Leak-proof workstation cover that doubles as endoscope transport bag to reduce biohazard risks
- Complies with guidelines for scope transport
- Drawstring ensures safety by keeping soiled scope securely enclosed during transfer
- Efficient padding to absorb moisture and protect scope during transport



Order Code	Description	Unit
BKTB	Single Use Scope Sheet and Transport Bag	50 / BOX

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
medicalcustomerservice@crkennedy.com.au | crkennedy.com.au

## Scope Transport Tray with Centre Reservoir


- Single-use rigid containment
- Bedside cleaning is performed within the Oasis Tray, keeping residual run-off contained
- Rigid construction protects delicate scopes, while the reversible lid clearly differentiates clean and soiled scopes
- Single-use alternative that protects both patients and staff from contamination and infection
- Oasis Trays are latex-free and biodegradable



Order Code	Description	Unit
OT1230	Oasis Scope Transport Tray with centre reservoir and reversible clean/soiled lid.	30/box

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
[medicalcustomerservice@crkennedy.com.au](mailto:medicalcustomerservice@crkennedy.com.au) | [crkennedy.com.au](http://crkennedy.com.au)



C.R. Kennedy  
MEDICAL DIVISION

ERCP





## FujiFilm® Medwork rotaCRUSH Lithotripter Extraction Basket

- Ready to use immediately without assembly
- Precise one-handed control of rotation
- Continuously adjustable basket size
- Proven Twist 'n' Catch basket that captures stones quickly and reliably
- Crown head design of the metal coil ensures precision force transmission to the stone and protects the basket from detachment
- The handle offers an integrated, additional power function for really tough cases
- Manufactured in Germany

Order Code	Description	Ø Basket	Working Length	Ø Working Channel	Unit
E502165	rotaCRUSH Lithotripter ERCP Basket, Integrated Catheter, Spiral Spring, Handle	35mm	200cm	3.8mm	Individual




Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
[medicalcustomerservice@crkennedy.com.au](mailto:medicalcustomerservice@crkennedy.com.au) | [crkennedy.com.au](http://crkennedy.com.au)



## FujiFilm® Medwork Twist 'n' Catch Stone Extraction Basket




- Maximum expansion force
- Precision control by rotary motion
- Infinite adjustment of basket size
- Proven Twist 'n' Catch basket captures stones quickly and reliably
- Spherical 130° rotation for short and long wire technique
- Manufactured in Germany

Order Code	Description	Ø Basket	Guide Wire	Working Length	Ø Working Channel	Unit
E502285	Twist 'n' Catch ECRP Basket for short and long wire technique	35mm	.035" ≥ 200cm	200cm	3.2mm	5/Box

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
[medicalcustomerservice@crkennedy.com.au](mailto:medicalcustomerservice@crkennedy.com.au) | [crkennedy.com.au](http://crkennedy.com.au)




C.R. Kennedy  
MEDICAL DIVISION

# Biopsy Forceps



## Designed to maximise sampling and streamline procedures

- Maximise your tissue sampling with the high volume cup
- Catheter markers assist with quick insertion and withdrawal
- Hinged jaws aid positioning at the sampling site
- Responsive closing force in tortuous anatomy



Catheter Markers



Serrated Cup

Order Code	Description	Cup Type	Ø Cup	Working Length	Ø Working Channel	Unit
WS-2423BUZ	Disposable Biopsy Forceps	Serrated	2.4mm	230cm	2.8mm	10/Box
WS-2423BNZ	Disposable Biopsy Forceps with needle	Serrated	2.4mm	230cm	2.8mm	10/Box
WS-2423BTZ	Disposable Biopsy Forceps	Oval	2.4mm	230cm	2.8mm	10/Box
WS-2423BWZ	Disposable Biopsy Forceps with needle	Oval	2.4mm	230cm	2.8mm	10/Box
WS-2423DG	Disposable Hot Biopsy Forceps	Oval	2.4mm	230cm	2.8mm	10/Box
WS-3023BTZ	Disposable Jumbo Biopsy Forceps	Oval	3.0mm	230cm	3.2mm	10/Box
WS-1816BT	Disposable Mini Biopsy Forceps	Oval	1.8mm	160cm	2.0mm	10/Box
WS-1810BT	Disposable Bronchoscope Biopsy Forceps	Oval	1.8mm	100cm	2.0mm	10/Box

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544  
[medicalcustomerservice@crkennedy.com.au](mailto:medicalcustomerservice@crkennedy.com.au) | [crkennedy.com.au](http://crkennedy.com.au)

Orders should be directed to C.R. Kennedy's Medical Division by email, fax, or phone.

Use the Product Order Code and include the quoted price and Quote Number (where available).

[medicalcustomerservice@crkennedy.com.au](mailto:medicalcustomerservice@crkennedy.com.au) | [www.crkennedy.com.au/medical](http://www.crkennedy.com.au/medical)

## SHIPMENT AND TERMS

Stock items are usually shipped within 48 hours of receiving an order.

Special orders for non-stock items will be subject to longer delivery times. There are no shipping charges for regular shipping. Urgent deliveries of stock items can be arranged by phone through your local C.R. Kennedy Medical Office. Urgent deliveries may be subject to additional delivery charges.

Payment terms are 30 days net from the end of the invoice month

## RETURNS POLICY

C.R. Kennedy will accept return of any endoscopy accessory stock item within 30 days of the original invoice date. Return authorization is required prior to return and all returned products must be in their original, undamaged and unopened packaging. Returned items are subject to a restocking charge and a credit will be applied to your account for future purchases. Non-stock items are not able to be returned under this returns policy

## C.R. KENNEDY MEDICAL OFFICE LOCATIONS

### Head Office

300 Lorimer St  
Port Melbourne VIC 3207  
Tel (03) 9823 1555  
Fax (03) 9827 7216

### New South Wales

Level 2, 15 Bourke Road  
Mascot NSW 2020  
Tel (02) 9552 8380  
Fax (02) 9552 8389

### Queensland

80 Kingsford-Smith Drive  
Albion QLD 4010  
Tel (07) 3632 6717  
Fax (07) 3632 6718

### South Australia

77 Fullarton Road  
Kent Town SA 5067  
Tel (08) 8410 1366  
Fax (08) 8410 1370

### Western Australia

17-19 Oxford Close  
West Leederville WA 6007  
Tel (08) 9381 7144  
Fax (08) 9381 7344

Contact C.R. Kennedy for more information:

VIC 03 9823 1812 | NSW 02 9552 8341 | QLD 07 3632 6717 | SA 08 8410 8105 | WA 08 9489 8544

[medicalcustomerservice@crkennedy.com.au](mailto:medicalcustomerservice@crkennedy.com.au) | [crkennedy.com.au](http://crkennedy.com.au)